

Product Name: Cytokeratin 7 Rabbit Monoclonal Antibody**Catalog #: AMRe86708**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.09mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:1000-1:5000,ICC/IF 1:100-1:200,FC 1:200-1:1000
Molecular Weight	Calculated MW:51 kDa; Observed MW:51 kDa

Antigen Information


Gene Name	Cytokeratin 7
Alternative Names	K7; Krt2-7; D15Wsu77e
Gene ID	110310
SwissProt ID	Q9DCV7
Immunogen	Recombinant protein of mouse Cytokeratin 7

Background

Blocks interferon-dependent interphase and stimulates DNA synthesis in cells.

Research Area

Image Data



Western blot analysis of extracts from MCF-7 cells using Cytokeratin 7 Rabbit Monoclonal Antibody at 1:1000.