

Product Name: LYPD3 Rabbit Monoclonal Antibody**Catalog #: AMRe86735**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:1000,IP 1:10-1:100
Molecular Weight	Calculated MW:36 kDa; Observed MW:75 kDa

Antigen Information


Gene Name	LYPD3
Alternative Names	C4.4A
Gene ID	27076
SwissProt ID	O95274
Immunogen	A synthetic peptide of human LYPD3

Background

Supports cell migration. May be involved in urothelial cell-matrix interactions. May be involved in tumor progression.

Research Area

Image Data



Western blot analysis of extracts from BxPC-3 cells using LYPD3 Rabbit Monoclonal Antibody at 1:1000.