
Product Name: ACTL8 Rabbit Monoclonal Antibody**Catalog #: AMRe86802**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,FC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:5000,FC 1:100-1:200,IP 1:50-1:100
Molecular Weight	Calculated MW:41 kDa; Observed MW:41 kDa


Antigen Information

Gene Name	ACTL8
Alternative Names	CT57
Gene ID	81569
SwissProt ID	Q9H568
Immunogen	A synthetic peptide of human ACTL8

Background

Research Area

Image Data



Western blot analysis of extracts from HepG2 cells using ACTL8 Rabbit Monoclonal Antibody at 1:1000.