
Product Name: RBM6 Rabbit Monoclonal Antibody**Catalog #: AMRe86819**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,ICC/IF 1:100-1:200
Molecular Weight	Calculated MW:129 kDa; Observed MW:129 kDa

Antigen Information


Gene Name	RBM6
Alternative Names	3G2; g16; DEF3; DEF-3; HLC-11; NY-LU-12
Gene ID	10180
SwissProt ID	P78332
Immunogen	A synthetic peptide of human RBM6

Background

Specifically binds poly(G) RNA homopolymers in vitro.

Research Area

Image Data



Western blot analysis of extracts from A549 cells using RBM6 Rabbit Monoclonal Antibody at 1:1000.