
Product Name: CD41 Rabbit Monoclonal Antibody**Catalog #: AMRe86831**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:500
Molecular Weight	Calculated MW:113 kDa; Observed MW:113 kDa

Antigen Information

Gene Name	CD41
Alternative Names	CD41; CD41B; GpIIb; AI172977; alphaIIb
Gene ID	16399
SwissProt ID	Q9QUM0
Immunogen	Recombinant protein of mouse CD41


Background

Integrin alpha-IIb/beta-3 is a receptor for fibronectin, fibrinogen, plasminogen, prothrombin, thrombospondin and vitronectin. It recognizes the sequence R-G-D in a wide array of ligands. It recognizes the sequence H-H-L-G-G-G-A-K-Q-A-G-D-V in

fibrinogen gamma chain. Following activation integrin alpha-IIb/beta-3 brings about platelet/platelet interaction through binding of soluble fibrinogen. This step leads to rapid platelet aggregation which physically plugs ruptured endothelial cell surface.

Research Area

Image Data



Western blot analysis of extracts from Mouse spleen tissue using CD41 Rabbit Monoclonal Antibody at 1:1000.