

Product Name: TRAFD1 Rabbit Monoclonal Antibody**Catalog #: AMRe86834**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:50-1:100
Molecular Weight	Calculated MW:65 kDa; Observed MW:65 kDa

Antigen Information

Gene Name	TRAFD1
Alternative Names	FLN29
Gene ID	10906
SwissProt ID	O14545
Immunogen	Recombinant protein of human TRAFD1


Background

The innate immune system confers host defense against viral and microbial infection, and TRAFD1 is a negative feedback regulator that controls excessive immune responses (Sanada et al., 2008 [PubMed 18849341]).[supplied by OMIM, Dec 2009]

Research Area

Image Data

Human spleen



Western blot analysis of extracts from Human spleen tissue using TRAFD1 Rabbit Monoclonal Antibody at 1:1000.