

---

**Product Name: CD133 Rabbit Monoclonal Antibody****Catalog #: AMRe86843**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,FC,IP
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:2000,IHC 1:200-1:1000,FC 1:50-1:200,IP 1:50-1:100
<b>Molecular Weight</b>	Calculated MW:97 kDa; Observed MW:133 kDa

**Antigen Information**

<b>Gene Name</b>	CD133
<b>Alternative Names</b>	RP41; AC133; CD133; MCDR2; STGD4; CORD12; PROML1; MSTP061
<b>Gene ID</b>	8842
<b>SwissProt ID</b>	O43490
<b>Immunogen</b>	Recombinant protein of human CD133


**Background**

This gene encodes a pentaspan transmembrane glycoprotein. The protein localizes to membrane protrusions and is often expressed on adult stem cells, where it is thought to function in maintaining stem cell properties by suppressing differentiation.

Mutations in this gene have been shown to result in retinitis pigmentosa and Stargardt disease. Expression of this gene is also associated with several types of cancer. This gene is expressed from at least five alternative promoters that are expressed in a tissue-dependent manner. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2009]

## Research Area

## Image Data



Western blot analysis of extracts from HT-29 cells using CD133 Rabbit Monoclonal Antibody at 1:1000.