

Product Name: HSPA12A Rabbit Monoclonal Antibody**Catalog #: AMRe86876**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:50-1:200,ICC/IF 1:100-1:200,IP 1:50-1:100
Molecular Weight	Calculated MW:75 kDa; Observed MW:75 kDa

Antigen Information


Gene Name	HSPA12A
Alternative Names	KIAA0417
Gene ID	259217
SwissProt ID	O43301
Immunogen	A synthetic peptide of human HSPA12A

Background

Heat shock 70 kDa protein 12A

Research Area

Image Data



Western blot analysis of extracts from Mouse spleen tissue using HSPA12A Rabbit Monoclonal Antibody at 1:1000.