

**Product Name: TAOK1 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe86894**



## Summary

<b>Production Name</b>	TAOK1 Rabbit Monoclonal Antibody
<b>Description</b>	Rabbit Monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.
<b>Purification</b>	Affinity Purification

## Immunogen

<b>Gene Name</b>	TAOK1
<b>Alternative Names</b>	PSK2; TAO1; KFC-B; MARKK; PSK-2; hKFC-B; hTAOK1; MAP3K16
<b>Gene ID</b>	57551
<b>SwissProt ID</b>	Q7L7X3.

## Application

<b>Dilution Ratio</b>	WB: 1:1000
<b>Molecular Weight</b>	Calculated MW:116 kDa; Observed MW:116 kDa

## Background


**Product Name: TAOK1 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe86894**




Serine/threonine-protein kinase involved in various processes such as p38/MAPK14 stress-activated MAPK cascade, DNA damage response and regulation of cytoskeleton stability

## Research Area

## Image Data



Western blot detection of TAOK1 in C6,Ramos using TAOK1 antibody(1:1000 diluted)



Western blot analysis of extracts from HepG2, RAW264.7 cells and Rat spleen tissue using AMRe86894 at 1:1000.

## Note

For research use only.