

Product Name: C2orf69 Rabbit Monoclonal Antibody**Catalog #: AMRe86916**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	Calculated MW:43 kDa; Observed MW:43 kDa

Antigen Information


Gene Name	C2orf69
Alternative Names	UPF0565 protein C2orf69
Gene ID	205327
SwissProt ID	Q8N8R5
Immunogen	Recombinant protein of human C2orf69

Background

UPF0565 protein C2orf69

Research Area

Image Data



Western blot analysis of extracts from U937 cells using C2orf69 Rabbit Monoclonal Antibody at 1:1000.