

Product Name: TMEM192 Rabbit Monoclonal Antibody
Catalog #: AMRe86957



Summary

Production Name	TMEM192 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, FC
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Immunogen

Gene Name	TMEM192
Alternative Names	Transmembrane protein 192
Gene ID	201931
SwissProt ID	Q8IY95.

Application

Dilution Ratio	WB: 1:1000-1:5000 ICC/IF: 1:100 FC: 1:100-1:200
Molecular Weight	Calculated MW:31 kDa; Observed MW:31 kDa

Background


Product Name: TMEM192 Rabbit Monoclonal Antibody
Catalog #: AMRe86957




Transmembrane protein 192

Research Area

Image Data



Western blot detection of TMEM192 in K562, HeLa cell lysates using TMEM192 antibody (1:1000 diluted).



Western blot analysis of extracts from SCC-9 cells using db14294 at 1:5000.

Note

For research use only.