

Product Name: TMEM161B Rabbit Monoclonal Antibody**Catalog #: AMRe86970**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:5000,IP 1:20-1:50
Molecular Weight	Calculated MW:56 kDa; Observed MW:56 kDa

Antigen Information


Gene Name	TMEM161B
Alternative Names	FLB3342; PRO1313
Gene ID	153396
SwissProt ID	Q8NDZ6
Immunogen	Recombinant protein of human TMEM161B

Background


Transmembrane protein 161B

Research Area

Image Data



Western blot detection of TMEM161B in K562,Hela cell lysates using TMEM161B antibody(1:1000 diluted).



Western blot analysis of extracts from T24 cells using AMRe86970 at 1:1000.