

Product Name: Flotillin 2 Rabbit Monoclonal Antibody
Catalog #: AMRe87010



Summary

Production Name	Flotillin 2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Immunogen

Gene Name	Flotillin 2
Alternative Names	ESA; ECS1; ESA1; ECS-1; M17S1
Gene ID	2319
SwissProt ID	Q14254.

Application

Dilution Ratio	WB: 1:1000-1:5000 IP: 1:20-1:50
Molecular Weight	Calculated MW:47 kDa; Observed MW:47 kDa

Background

Caveolae are small domains on the inner cell membrane involved in vesicular trafficking and signal transduction. This gene

Product Name: Flotillin 2 Rabbit Monoclonal Antibody
Catalog #: AMRe87010



encodes a caveolae-associated, integral membrane protein, which is thought to function in neuronal signaling. [provided by RefSeq, Jul 2008]

Research Area

Image Data



Western blot analysis of extracts from HeLa cells using AMRe87010 at 1:1000.

Note

For research use only.