

Product Name: POC5 Rabbit Monoclonal Antibody**Catalog #: AMRe87034**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	Calculated MW:63 kDa; Observed MW:70 kDa

Antigen Information


Gene Name	POC5
Alternative Names	C5orf37
Gene ID	134359
SwissProt ID	Q8NA72
Immunogen	Recombinant protein of human POC5

Background


Essential for the assembly of the distal half of centrioles, required for centriole elongation.

Research Area

Image Data



Western blot analysis of extracts from HeLa cells using POC5 Rabbit Monoclonal Antibody at 1:1000.



Western blot analysis of extracts from HCT116 cells using AMRe87034 at 1:1000.