

---

**Product Name: Pirh2 Rabbit Monoclonal Antibody****Catalog #: AMRe87039**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ICC/IF
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:2000,IHC 1:50-1:100,ICC/IF 1:100-1:200
<b>Molecular Weight</b>	Calculated MW:30 kDa; Observed MW:30 kDa

**Antigen Information**

<b>Gene Name</b>	Pirh2
<b>Alternative Names</b>	ZCHY; ARNIP; CHIMP; PIRH2; RNF199; ZNF363; PRO1996
<b>Gene ID</b>	25898
<b>SwissProt ID</b>	Q96PM5
<b>Immunogen</b>	Recombinant protein of human Pirh2


**Background**

The protein encoded by this gene has ubiquitin ligase activity. It mediates E3-dependent ubiquitination and proteasomal degradation of target proteins, including tumor protein 53, histone deacetylase 1, and cyclin-dependent kinase inhibitor 1B,

thus regulating their levels and cell cycle progression. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jun 2013]

## Research Area

## Image Data



Western blot analysis of extracts from HeLa cells using Pirh2 Rabbit Monoclonal Antibody at 1:1000.