

Product Name: CD31 Rabbit Monoclonal Antibody**Catalog #: AMRe87064**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IP 1:20-1:50
Molecular Weight	Calculated MW:81 kDa; Observed MW:130 kDa

Antigen Information

Gene Name	CD31
Alternative Names	Cd31; Pecam; C85791; PECAM-1
Gene ID	18613
SwissProt ID	Q08481
Immunogen	Recombinant protein of mouse CD31

Background


Cell adhesion molecule which is required for leukocyte transendothelial migration (TEM) under most inflammatory conditions (By similarity). Tyr-679 plays a critical role in TEM and is required for efficient trafficking of PECAM1 to and from the lateral

border recycling compartment (LBRC) and is also essential for the LBRC membrane to be targeted around migrating leukocytes (By similarity).

Research Area

Image Data

Mouse lung



kDa

250 -

150 -

100 - }
75 -

50 -

37 -

25 -

20 -

15 -

10 -

Western blot analysis of extracts from Mouse lung tissue using CD31 Rabbit Monoclonal Antibody at 1:1000.