
Product Name: 53BP1 Rabbit Monoclonal Antibody**Catalog #: AMRe87074**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:5000-1:50000,IHC 1:50-1:100,ICC/IF 1:100-1:200
Molecular Weight	Calculated MW:214 kDa; Observed MW:450 kDa

Antigen Information

Gene Name	53BP1
Alternative Names	TP53; p202; 53BP1; TDRD30; p53BP1
Gene ID	7158
SwissProt ID	Q12888
Immunogen	A synthetic peptide of human 53BP1


Background

This gene encodes a protein that functions in the DNA double-strand break repair pathway choice, promoting non-homologous end joining (NHEJ) pathways, and limiting homologous recombination. This protein plays multiple roles in the

DNA damage response, including promoting checkpoint signaling following DNA damage, acting as a scaffold for recruitment of DNA damage response proteins to damaged chromatin, and promoting NHEJ pathways by limiting end resection following a double-strand break. These roles are also important during V(D)J recombination, class switch recombination and at unprotected telomeres. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2017]

Research Area

Image Data



Western blot analysis of extracts from HeLa cells using 53BP1 Rabbit Monoclonal Antibody at 1:1000.