
Product Name: Zinc Alpha 2 Glycoprotein Rabbit Monoclonal Antibody**Catalog #: AMRe87145**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:100-1:200,IP 1:20-1:50
Molecular Weight	Calculated MW:34 kDa; Observed MW:41 kDa

Antigen Information


Gene Name	Zinc Alpha 2 Glycoprotein
Alternative Names	ZAG; ZA2G
Gene ID	563
SwissProt ID	P25311
Immunogen	A synthetic peptide of human Zinc Alpha 2 Glycoprotein

Background

Stimulates lipid degradation in adipocytes and causes the extensive fat losses associated with some advanced cancers

Research Area

Image Data



Western blot analysis of extracts from HepG2 cells and Human liver tissue using Zinc Alpha 2 Glycoprotein Rabbit Monoclonal Antibody at 1:1000.