

Product Name: p53 DINP1/TP53INP1 Rabbit Monoclonal Antibody**Catalog #: AMRe87161**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:100-1:200,ICC/IF 1:100-1:200
Molecular Weight	Calculated MW:27 kDa; Observed MW:27 kDa

Antigen Information


Gene Name	p53 DINP1/TP53INP1
Alternative Names	SIP; Teap; p53DINP1; TP53DINP1; TP53INP1A; TP53INP1B
Gene ID	94241
SwissProt ID	Q96A56
Immunogen	A synthetic peptide of human p53 DINP1/TP53INP1

Background

Antiproliferative and proapoptotic protein involved in cell stress response which acts as a dual regulator of transcription and autophagy

Research Area

Image Data



Western blot analysis of extracts from 3T3 cells using p53 DINP1/TP53INP1 Rabbit Monoclonal Antibody at 1:1000.