
Product Name: SerpinB2/PAI-2 Rabbit Monoclonal Antibody**Catalog #: AMRe87187**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IP 1:20-1:50
Molecular Weight	Calculated MW:47 kDa; Observed MW:47 kDa

Antigen Information


Gene Name	SerpinB2/PAI-2
Alternative Names	PAI; PAI2; PAI-2; PLANH2; HsT1201
Gene ID	5055
SwissProt ID	P05120
Immunogen	Recombinant protein of human SerpinB2/PAI-2

Background

Inhibits urokinase-type plasminogen activator

Research Area

Image Data



Western blot analysis of extracts from U937 cells using SerpinB2/PAI-2 Rabbit Monoclonal Antibody at 1:1000.