
Product Name: MIOX Rabbit Monoclonal Antibody**Catalog #: AMRe87237**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IP 1:100-1:200
Molecular Weight	Calculated MW:33 kDa; Observed MW:33 kDa

Antigen Information

Gene Name	MIOX
Alternative Names	ALDRL6
Gene ID	55586
SwissProt ID	Q9UGB7
Immunogen	A synthetic peptide of human MIOX

Background

Fe cationSearch proteins in UniProtKB for this molecule

Research Area

Image Data

Human fetal kidney

kDa

250 -

150 -

100 -

75 -

50 -

37 -

25 -

20 -

15 -

10 -



Western blot analysis of extracts from Human fetal kidney tissue using MIOX Rabbit Monoclonal Antibody at 1:1000.