

Product Name: TPPP Rabbit Monoclonal Antibody**Catalog #: AMRe87270**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:100-1:200
Molecular Weight	Calculated MW:24 kDa; Observed MW:24 kDa

Antigen Information


Gene Name	TPPP
Alternative Names	p24; p25; TPPP1; TPPP/p25; p25alpha
Gene ID	11076
SwissProt ID	O94811
Immunogen	A synthetic peptide of human TPPP

Background

Regulator of microtubule dynamics that plays a key role in myelination by promoting elongation of the myelin sheath (PubMed:31522887)

Research Area

Image Data



Western blot analysis of extracts from SH-SY5Y cells using TPPP Rabbit Monoclonal Antibody at 1:1000.