

Product Name: SerpinB2/PAI-2 Rabbit Monoclonal Antibody**Catalog #: AMRe87281**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.15mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:500-1:2000
Molecular Weight	Calculated MW:47 kDa; Observed MW:47 kDa

Antigen Information


Gene Name	SerpinB2/PAI-2
Alternative Names	PAI; PAI2; PAI-2; PLANH2; HsT1201
Gene ID	5055
SwissProt ID	P05120
Immunogen	Recombinant protein of human SerpinB2/PAI-2

Background

Inhibits urokinase-type plasminogen activator

Research Area

Image Data



Western blot analysis of extracts from T24 cells using SerpinB2/PAI-2 Rabbit Monoclonal Antibody at 1:1000.