
Product Name: PRA1 Rabbit Monoclonal Antibody**Catalog #: AMRe87309**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:10000
Molecular Weight	Calculated MW:21 kDa; Observed MW:21 kDa

Antigen Information


Gene Name	PRA1
Alternative Names	PRA1; YIP3; PRAF1
Gene ID	10567
SwissProt ID	Q9UI14
Immunogen	A synthetic peptide of human PRA1

Background

General Rab protein regulator required for vesicle formation from the Golgi complex

Research Area

Image Data



Western blot analysis of extracts from HeLa cells using PRA1 Rabbit Monoclonal Antibody at 1:1000.