

Product Name: HDLBP Rabbit Monoclonal Antibody**Catalog #: AMRe87335**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:2000-1:20000,IHC 1:100-1:200
Molecular Weight	Calculated MW:142 kDa; Observed MW:142 kDa

Antigen Information


Gene Name	HDLBP
Alternative Names	HBP; VGL; PRO2900
Gene ID	3069
SwissProt ID	Q00341
Immunogen	A synthetic peptide of human HDLBP

Background


The protein encoded by this gene binds high density lipoprotein (HDL) and may function to regulate excess cholesterol levels in cells. The encoded protein also binds RNA and can induce heterochromatin formation. [provided by RefSeq, Mar 2016]

Research Area

Image Data



Western blot detection of HDLBP in HeLa, A549, HL-60, U251, U87-MG cell lysates using HDLBP antibody (1:1000 diluted).



Western blot analysis of extracts from HeLa, C2C12 cells and Mouse muscle, Rat heart tissue using AMRe87335 at 1:5000.