
Product Name: Phospho-p130 (Thr986) Rabbit Monoclonal Antibody**Catalog #: AMRe87345**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,FC,IP
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,FC 1:100-1:200,IP 1:50-1:200
Molecular Weight	Calculated MW:128 kDa; Observed MW:130 kDa

Antigen Information

Gene Name	Phospho-p130
Alternative Names	Rb2; P130
Gene ID	5934
SwissProt ID	Q08999
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Thr986 of human p130


Background

Key regulator of entry into cell division. Directly involved in heterochromatin formation by maintaining overall chromatin


structure and, in particular, that of constitutive heterochromatin by stabilizing histone methylation. Recruits and targets histone methyltransferases KMT5B and KMT5C, leading to epigenetic transcriptional repression. Controls histone H4 'Lys-20' trimethylation. Probably acts as a transcription repressor by recruiting chromatin-modifying enzymes to promoters. Potent inhibitor of E2F-mediated trans-activation, associates preferentially with E2F5. Binds to cyclins A and E. Binds to and may be involved in the transforming capacity of the adenovirus E1A protein. May act as a tumor suppressor.

Research Area

Image Data



Western blot detection of Phospho-Rb2 p130 (Thr986) in K562 cell lysates using Phospho-Rb2 p130 (Thr986) antibody (1:1000 diluted).



Western blot analysis of extracts from C2C12 cells and Mouse muscle tissue using AMRe87345 at 1:1000.