
Product Name: UCH37 Rabbit Monoclonal Antibody**Catalog #: AMRe87349**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:5000,IHC 1:200-1:1000
Molecular Weight	Calculated MW:38 kDa; Observed MW:38 kDa

Antigen Information


Gene Name	UCH37
Alternative Names	UCH37; CGI-70; INO80R; UCH-L5
Gene ID	51377
SwissProt ID	Q9Y5K5
Immunogen	A synthetic peptide of human UCH37

Background

Protease that specifically cleaves 'Lys-48'-linked polyubiquitin chains

Research Area

Image Data



Western blot analysis of extracts from Human fetal kidney tissue using UCH37 Rabbit Monoclonal Antibody at 1:1000.