

Product Name: CD8B Rabbit Monoclonal Antibody**Catalog #: AMRe87402**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:500,IP 1:20-1:50
Molecular Weight	Calculated MW:24 kDa; Observed MW:27-32 kDa

Antigen Information


Gene Name	CD8B
Alternative Names	Cd8b; Ly-3; Ly-C; Lyt-3
Gene ID	12526
SwissProt ID	P10300
Immunogen	Recombinant protein of mouse CD8B

Background

Integral membrane glycoprotein that plays an essential role in the immune response and serves multiple functions in responses against both external and internal offenses

Research Area

Image Data



Western blot analysis of extracts from Mouse thymus tissue using CD8B Rabbit Monoclonal Antibody at 1:1000.