

Product Name: IGFBP5 Rabbit Monoclonal Antibody**Catalog #: AMRe87404**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IP 1:20-1:50
Molecular Weight	Calculated MW:30 kDa; Observed MW:30 kDa

Antigen Information


Gene Name	IGFBP5
Alternative Names	IGFBP-5; AI256729; AW208790; IGFBP-5P
Gene ID	16011
SwissProt ID	P24594
Immunogen	Recombinant protein of mouse IGFBP5

Background

IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture

Research Area

Image Data



Western blot analysis of extracts from Rat heart tissue using IGFBP5 Rabbit Monoclonal Antibody at 1:1000.