

Product Name: uPA Receptor/U-PAR Rabbit Monoclonal Antibody**Catalog #: AMRe87420**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,FC,IP
Reactivity	Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,FC 1:200-1:500,IP 1:20-1:50
Molecular Weight	Calculated MW:35 kDa; Observed MW:35-65 kDa

Antigen Information


Gene Name	uPA Receptor/U-PAR
Alternative Names	Cd87; uPAR; u-PAR
Gene ID	18793
SwissProt ID	P35456
Immunogen	Recombinant protein of mouse uPA Receptor/U-PAR

Background

Acts as a receptor for urokinase plasminogen activator

Research Area

Image Data



Western blot analysis of extracts from Mouse lung tissue using uPA Receptor/U-PAR Rabbit Monoclonal Antibody at 1:1000.