
Product Name: SLAM/CD150 Rabbit Monoclonal Antibody**Catalog #: AMRe87421**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:500
Molecular Weight	Calculated MW:38 kDa; Observed MW:70-85 kDa

Antigen Information


Gene Name	SLAM/CD150
Alternative Names	Slam; CD150; IPO-3; CDw150; ESTM51; AA177906; 4933415F16
Gene ID	27218
SwissProt ID	Q9QUM4
Immunogen	Recombinant protein of mouse SLAM/CD150

Background

Self-ligand receptor of the signaling lymphocytic activation molecule (SLAM) family

Research Area

Image Data



Western blot analysis of extracts from Rat thymus tissue using SLAM/CD150 Rabbit Monoclonal Antibody at 1:1000.