

Product Name: CD134 Rabbit Monoclonal Antibody**Catalog #: AMRe87534**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,FC,IP
Reactivity	Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:1000,FC 1:200-1:500,IP 1:20-1:50
Molecular Weight	Calculated MW:31 kDa; Observed MW:45 kDa

Antigen Information


Gene Name	CD134
Alternative Names	Ox40; ACT35; CD134; Ly-70; Txgp1; TXGP1L
Gene ID	22163
SwissProt ID	P47741
Immunogen	Recombinant protein of mouse CD134

Background

Receptor for TNFSF4/OX40L/GP34. Is a costimulatory molecule implicated in long-term T-cell immunity (By similarity).

Research Area

Image Data



Western blot detection of CD134 in mouse spleen cell lysates using CD134 antibody(1:1000 diluted).