
Product Name: FOXL2 Rabbit Monoclonal Antibody**Catalog #: AMRe87538**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:2000,ICC/IF 1:20-1:50,IP 1:20-1:50
Molecular Weight	Calculated MW:39 kDa; Observed MW:49 kDa

Antigen Information

Gene Name	FOXL2
Alternative Names	BPES; PFRK; POF3; BPES1; PINTO
Gene ID	668
SwissProt ID	P58012
Immunogen	A synthetic peptide of human FOXL2


Background

This gene encodes a forkhead transcription factor. The protein contains a fork-head DNA-binding domain and may play a role in ovarian development and function. Expansion of a polyalanine repeat region and other mutations in this gene are a cause of

blepharophimosis syndrome and premature ovarian failure 3. [provided by RefSeq, Jul 2016]

Research Area

Image Data



Western blot analysis of extracts from HeLa cells using FOXL2 Rabbit Monoclonal Antibody at 1:1000.