

Product Name: Angiotensinogen Rabbit Monoclonal Antibody**Catalog #: AMRe87613**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:200-1:2000,IP 1:20-1:50
Molecular Weight	Calculated MW:52 kDa; Observed MW:52 kDa

Antigen Information


Gene Name	Angiotensinogen
Alternative Names	AGT;Angl; AngII; Aogen; AI265500; Serpina8
Gene ID	11606
SwissProt ID	P11859
Immunogen	Recombinant protein of mouse Angiotensinogen

Background

Essential component of the renin-angiotensin system (RAS), a potent regulator of blood pressure, body fluid and electrolyte homeostasis.

Research Area

Image Data



Western blot detection of Angiotensinogen in mouse liver cell lysates using Angiotensinogen antibody(1:1000 diluted).