

---

**Product Name: NMT55/p54nrb Rabbit Monoclonal Antibody****Catalog #: AMRe87627**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ICC/IF
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:1000-1:5000,IHC 1:200-1:500,ICC/IF 1:500-1:2000
<b>Molecular Weight</b>	Calculated MW:54 kDa; Observed MW:54 kDa

**Antigen Information**

<b>Gene Name</b>	NMT55/p54nrb
<b>Alternative Names</b>	P54; NMT55; NRB54; MRXS34; P54NRB; PPP1R114
<b>Gene ID</b>	4841
<b>SwissProt ID</b>	Q15233
<b>Immunogen</b>	A synthetic peptide of human NMT55/p54nrb


**Background**

This gene encodes an RNA-binding protein which plays various roles in the nucleus, including transcriptional regulation and RNA splicing. A rearrangement between this gene and the transcription factor E3 gene has been observed in papillary renal cell

carcinoma. Alternatively spliced transcript variants have been described. Pseudogenes exist on Chromosomes 2 and 16.  
[provided by RefSeq, Feb 2009]

## Research Area

## Image Data



Western blot analysis of extracts from RAW264.7, Jurkat cells and Rat spleen tissue using NMT55/p54nrb Rabbit Monoclonal Antibody at 1:2000.