

Product Name: RABL3 Rabbit Monoclonal Antibody**Catalog #: AMRe87779**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:50-1:100,ICC/IF 1:100-1:200,FC 1:50-1:200
Molecular Weight	Calculated MW:26 kDa; Observed MW:26 kDa

Antigen Information


Gene Name	RABL3
Alternative Names	Rab-like protein 3
Gene ID	285282
SwissProt ID	Q5HYI8
Immunogen	A synthetic peptide of human RABL3

Background


RAB, member of RAS oncogene family like 3

Research Area

Image Data



Western blot analysis of extracts from 293 cells using RABL3 Rabbit Monoclonal Antibody at 1:1000.



Western blot analysis of extracts from HeLa, HepG2, HCT116 cells using AMRe87779 at 1:1000.