

Product Name: Ly6A/E Rabbit Monoclonal Antibody
Catalog #: AMRe87810



Summary

Production Name	Ly6A/E Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, FC
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Immunogen

Gene Name	Ly6A/E
Alternative Names	TAP; Sca1; Sca-1; Ly-6A.2; Ly-6A/E; Ly-6E.1
Gene ID	110454
SwissProt ID	P05533.

Application

Dilution Ratio	WB: 1:1000 FC: 1:200-1:1000
Molecular Weight	Calculated MW:14 kDa; Observed MW:18 kDa

Background


T-cell activation

Product Name: Ly6A/E Rabbit Monoclonal Antibody
Catalog #: AMRe87810



Research Area

Image Data



Western blot analysis of extracts from Mouse kidney tissue using AMRe87810 at 1:1000.

Note

For research use only.