

Product Name: CRMP4 Rabbit Monoclonal Antibody
Catalog #: AMRe87818



Summary

Production Name	CRMP4 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.
Purification	Affinity Purification

Immunogen

Gene Name	CRMP4
Alternative Names	DRP3; ULIP; CRMP4; DRP-3; LCRMP; CRMP-4; ULIP-1
Gene ID	1809
SwissProt ID	Q14195.

Application

Dilution Ratio	WB: 1:1000-1:5000
Molecular Weight	Calculated MW:62 kDa; Observed MW:62 kDa

Background

Necessary for signaling by class 3 semaphorins and subsequent remodeling of the cytoskeleton


Product Name: CRMP4 Rabbit Monoclonal Antibody
Catalog #: AMRe87818



Research Area

Image Data

Human hippocampus



Western blot analysis of extracts from Human hippocampus tissue using AMRe87818 at 1:1000.

Note

For research use only.