

**Product Name: CHIKV Spike glycoprotein E2
Recombinant Antibody (CHK265)
Catalog #: AMRe60084**

Summary

Production Name	CHIKV Spike glycoprotein E2 Recombinant Antibody (CHK265)
Description	Recombinant Monoclonal antibody
Host	Mammalian cells
Application	ELISA, FC
Reactivity	Chikungunya virus (strain S27-African prototype) (CHIKV)

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	human IgG1, kappa
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	0.01M PBS, pH 7.4.
Purification	Affinity Purified

Immunogen

Gene Name	E2
Alternative Names	Structural polyprotein, Envelope glycoprotein E2, Spike glycoprotein E2
SwissProt ID	Q8JUX5.

Application

Dilution Ratio	ELISA 1:1000-2000 FC 1:50-100
-----------------------	-------------------------------

Background

Chikungunya virus (CHIKV) is an arthropod-borne alphavirus that causes Chikungunya fever, a re-emerging disease characterized by high fever, rash, and arthralgia. The E2 envelope glycoprotein is a structural protein that forms heterodimers with the E2 glycoprotein and is displayed on the viral surface as part of trimeric spike complexes. E2 is primarily



**Product Name: CHIKV Spike glycoprotein E2
Recombinant Antibody (CHK265)
Catalog #: AMRe60084**

responsible for receptor binding and host cell attachment, playing a key role in viral entry.

Note

For research use only.