

Product Name: Anti-Langya virus/LayV F/Fusion Protein Antibody (SAA2264)

Catalog #: AMM60148

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	ELISA, WB
Reactivity	Mojiang virus
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.

Application

Dilution Ratio ELISA 1:5000-1:20000 WB 1:500-1:2000

Molecular Weight

Antigen Information

Gene Name F/Fusion Protein
Alternative Names F/Fusion Protein
Gene ID
SwissProt ID
Immunogen

Background

Nipah virus (NiV), a highly pathogenic zoonotic paramyxovirus, mediates host cell entry through the coordinated action of two envelope glycoproteins: the attachment glycoprotein (G), which engages the host cell receptors ephrin-B2 and ephrin-B3, and the fusion glycoprotein (F), which drives membrane fusion following receptor binding. These viral surface proteins, along with

the viral RNA-dependent RNA polymerase (RdRp) complex comprising the nucleoprotein (N), phosphoprotein (P), and large polymerase protein (L), constitute the primary molecular targets for therapeutic intervention and vaccine development. The conserved receptor-binding interface of the G protein and the fusion machinery of the F protein represent particularly attractive targets for neutralizing antibodies and entry inhibitors, while the viral replication machinery offers potential targets for broad-spectrum antiviral compounds targeting paramyxovirus transcription and replication.

Research Area

Microbiology&Infectious Disease